

Lining

Inside Lining of PPGL / PFRP / S.STEEL (316 / 316L / 304) ORAC Cement sheet dully epoxy coated as acid

Work Space

Table work space provided with Acid / Alkali resistance Tiles / Jet Black Granite / Laminate / S. Steel / PP FRP

Services

Ceramic / SS / PVC Sink, Water Line with Swan neck tap, gas / vacuum line etc.....

Under Bench

In model offered with Under bench is having provision for storage and LPG Cylinder storing with middle shelf

Facia

Electrical points of 5 / 15 Amp. Socket, Sandwich for illumination. Remote Control Knobs for Water / LPG /

Exhaust

Exhaust system comprising of centrifugal blower made of PPGL with FRP lining with dynamically balanced impeller

Duct & Fittings

Generally 3MM thick plain or 6 mm thick PPGL Lined ducting of suitable size for inlet & outlet provided with valves

(*Optional).

Due to continuous R & D and specifications may change. Pictures shown with Optional accessories.

Table Top Model

900 X 750 X 1800mm

1200 X 750 X 1800mm

1500 X 750 X 1800mm

Model with under Bench

Storing

900 X 750 X 2700mm

1200 X 750 X 2700mm

1500 X 750 X 2700mm

Customized size on request