

#### QUALITY ACCESSORIES :

Hardware like hinges, handles etc on panels are from best indigenous / foreign suppliers.

Optional accessories like : Digital Temperature Controller / Indicator, Audio Temperature Alarm, Strip Curtains, Inner Storage Racks etc.....are available on request at extra.

#### TEMP./HUMIDITY RANGES :

Temp. Ranges : As per Customer requirement (-20 C, -10C, 0C, + 50C & Above ambient. )

Humidity Ranges : 40% to 95% RH (optional for walk-in-humidity chambers)

#### DATA RECORDING FACILITIES

- 1) Circular (Daily/Weekly) or Strip Chart (Continues) Recorder.
- 2) Printer Interface compitable to 80 column DOT Matrix printer (EPSON/WIPTRO)
- 3) PC Interface through RS 232 - C/RS 485 for multiple networking of controllers.
- 4) Window based Software for reports in Tubular / Tabular/Graphical form.
- 5) Validation : HMG offers validation protocols for DQ, IQ, & PQ documentation.